




Rosenmund Filter[®]

Side Discharge

RSD 8910 US



ROSENMUND

Rosenmund Filter® Side Discharge


Benefits

The Rosenmund Filter® is installed for solid-liquid separation of suspensions and for washing and isolation of solids in the chemical, fine chemical and pharmaceutical processing industries. It can be used either for a dedicated product or as a multi-purpose unit.


The side discharge filter is available as a filter-dryer, see bulletin RFT-S.

The side discharge filter is ideal for the following conditions:


- frequent changes of product
- processes calling for the highest standards of particle-freedom and sterility
- limited clear height of building




Charging and filtration:
vacuum or pressure



Washing:
re-slurry or displacement



Smoothing
and compressing




Discharging

Features

The Rosenmund Filter® with side discharge is equipped with a specially contoured beam-agitator and a patented metal sealing side discharge valve.

- The side discharge valve permits pressure tight shut-off after every product discharge without previous cleaning of the sealing surfaces, a distinct advantage in automated processing and contamination-free production.
- The special impeller (two- or three-armed) facilitates uniform agitation and discharge. The three-armed version has in addition a strong centering effect on the agitation system for improved stability. Both designs ensure a quiet and easily installed and maintained operation.

- ASME code vessels, 50 psig or 90 psig to full vacuum, in all grades of stainless steel, hastelloy and other alloys, are locally manufactured in the United States.
- Hydraulic power pack, electric or hydraulic motors and control system.
- Fully trained, locally based, Rosenmund service engineers are available for contract maintenance.
- Complete turnkey process systems, skids, and custom designs are available from Rosenmund Engineering Services.
- Laboratory and pilot test units are also available.
- Complete stock of spare parts, locally available.



Schematic diagram

Filter- Size m ²	Filtra- tion Area sq. ft.	Solids Vol. cu. ft.	Max. Liquid Vol. gal.	D in.	H in.	S in.
1	11	17	290	47	131	20
1,6	17	28	450	59	131	20
2,5	27	44	750	71	134	20
3	32	53	900	79	150	20
4	43	71	1400	94	156	20
5	54	88	1700	102	157	20
6	65	107	2000	110	159	20
8	86	170	3000	130	183	24
10	108	214	4100	146	187	24

ROSENMUND INC.

2969 Interstate Street

P.O. Box 668625

Charlotte NC 28266-8625

Telephone (704) 398-1111

Fax (704) 398 1089


ROSENMUND INC.

Filtration, Drying, Vacuum Equipment & Engineered Systems


INSTALLATION & EQUIPMENT HIGHLIGHTS

3M² Rosenmund Filter[®] Side Discharge

For details, please see reverse side.



Installed at Syntex Chemicals Inc., Boulder, Colorado



Size	3M ²
Type	Side Discharge
Filtration Area	32 ft. ²
Working Slurry Volume	725 gallons
Solids Volume	56 ft. ³
Design Pressure	50 psig/full vacuum
Design Temperature	302°F
Material of Construction	316L Stainless Steel
	ASME Code
Agitator	Three-arm, S-Blade
Solids Discharge	Side, hydraulic valve
Side Discharge Valve	Rosenmund patented metal to metal sealing
Filter Plate	Simple heated
	90 psig, 302°F
Jacket	Half-pipe coil
	125 psig, 302°F
Temperature Measurement	Thermowell
Dust Filter	Blowback, bag-type
Motors	30 HP, hydraulic, agitator rotation
	5 HP, hydraulic, all raising/lowering operations and valve open/close
Weight (empty)	20,000 pounds
Installation Date	June, 1989

Application Details

Syntex Chemicals Inc., located in Boulder, Colorado, has installed this filter as part of its new multi-purpose manufacturing facility. The Rosenmund filter is used to produce pharmaceutical intermediates and final products. Filtration, displacement washing, re-slurry, drying, and automatic discharge are performed in the unit. Drying is accomplished by hot gas recirculation using an APOVAC® Vacuum/Gas Compressor System and by pulling vacuum either through the dust filter or the filtrate outlets. The simple heated filterplate and full jacket provides heat energy into the cake. The metal to metal side discharge valve allows continuous sealing without the necessity for between batch cleaning. To complete the installation, Rosenmund supplied the heat transfer module (HTM) and fully automated controls using an Allen Bradley PLC.

ROSENMUND

Rosenmund Inc.
P.O. Box 668625
Charlotte, N.C. 28266-8625
704-398-1111 (Tel) 704-398-1089 (Fax)

Rosenmund Inc. serves the pharmaceutical, chemical, specialty chemical, and process industries with a full range of filters, filter-dryers, and Apovac® vacuum equipment. Engineered systems include solvent recovery, sterile and CIP designs, recirculation drying, heat transfer, reactor to filter-dryer packages, and customized applications.

Rosenmund Filter-Dryer®

Center Discharge – Heated Paddle System

RFT-P 9004 US



ROSENMUND

Rosenmund Filter-Dryer® Center Discharge - Heated Paddle System

Applications


Rosenmund Filter-Dryers® can be employed wherever a drying stage follows filtration. According to application requirements, either contact drying (under vacuum) or convection drying (e.g. with closed nitrogen circuit) can be selected as the drying method.

The Rosenmund Filter® is ideally suited for drying. To provide maximum heat transfer into the product, three separate components of the filter are heated: vessel (wall, base, head), filter plate, and agitator system.

- The entire vessel can be jacketed with half-pipe coil for heating fluid or steam circulation.

- The filter plate is directly heated with steam or heating fluid.
- The agitator is heated by circulating hot fluid through aerofoil-designed paddles.

Complete turnkey process systems, skid packages (including heat transfer and vacuum/gas recirculation systems), and custom designs are available from Rosenmund Engineering Services. Pilot test equipment, for on-site testwork, provide detail characteristics and scale-up data for your process.



Benefits

- Filtration and drying in a single vessel
- Lower capital and maintenance costs
- Less plant space required
- Lower total installed cost
- Improved product quality with less product handling
- Single source responsibility

Key

- 1 Rosenmund Filter-Dryer® with connection for heating/cooling media, nozzles for slurry (a), wash-liquor (b), pressure gas (c), filtrate (d), product discharge (e), gas circulation (f)
- 2 Dust filter with heating jacket and nitrogen blowback system (g)
- 3 Pre-condenser
- 4 Ring liquid compressor/Vacuum pump
- 5 Ring liquid cooler
- 6 Ring liquid tank with condensate or overflow (h)
- 7 Exhaust gas /condenser with exhaust (i) and gas circulation connection (f)
- 8 Demister
- 9 Gas heater (for gas circulation)


Process Designs

Pressure filtration: The compressor can generate the required pressure differential for filtration.

Dry blowing: After filtration, gas is blown through the cake from top to bottom (in a closed gas loop) to remove additional moisture prior to drying.

Convection drying: The cake is blown through from above with warmed gas, resulting in convective drying of the cake.

Vacuum drying: Solvent vapors pass through the dust filter and condenser into the vacuum system. Alternatively, vacuum is pulled through the filtrate outlets. Non-condensable gases are vented to the atmosphere via a vent condenser. The gas is solvent-free to meet applicable volatile organic compound emission regulations.



Filter Size sq.ft	Solids Vol. cu.ft.	Max. Liquid Vol. gal.	D in	H in	L in	S in	A in	B in (approx.)
32	53	1200	79	222	41	20	120	90
43	71	1800	95	228	38	20	120	90
54	88	2200	102	230	38	20	120	90
65	106	2600	110	232	37	20	120	90
86	141	4000	130	258	36	20	132	102
108	176	4200	146	262	35	20	132	102
129	212	6000	157	266	35	20	132	102
172	282	10000	173	270	35	20	132	102

ROSENMUND INC.

2969 Interstate Street P.O. Box 668625 Charlotte NC 28266-8625 Telephone (704) 398-1111 Fax (704) 398 1089


ROSENMUND INC.


Filtration, Drying, Vacuum Equipment & Engineered Systems

INSTALLATION & EQUIPMENT HIGHLIGHTS

4M² Rosenmund Filter-Dryer[®]
Center Discharge-heated
paddle system

For details, please see reverse side.





Size	4M ²
Type	Paddle Filter-Dryer
Filtration Area	43 ft. ²
Working Slurry Volume	1200 gallons
Solids Volume	71 ft. ³
Design Pressure	90 psig/full vacuum
Design Temperature	-20° to 370°F
Material of Construction	316L stainless steel ASME Code
Agitator	Three-arm paddle, heated
Solids Discharge	Center
Filter Plate	Simple heated 150 psig, 370°F
Jacket	Half-pipe coil 150 psig, 370°F
Temperature Measurement	Thermowell
Dust Filter	20 in blanked nozzle for future installation
Motors	30 HP, hydraulic, agitator rotation 3 HP, hydraulic, all raising/lowering operations
Weight (empty)	24,000 pounds
Installation Date	July, 1990

Application Details

This filter-dryer replaced existing centrifuges and multi-leaf filters and, in the future, will eliminate existing tumble dryers. The filter-dryer is being employed for three applications: two are for catalyst cleanup and disposal, where the filtrate (liquid) is the product and one where the dried cake is the product. The unit is equipped with a slurry outlet valve at the filter plate level for slurry/liquid discharge. A Rosenmund supplied PLC control system operates the filter-dryer. Rosenmund also supplied a Kammerer Knife Gate Valve to provide for a pressure/vacuum tight seal at the discharge tube. Installation benefits include flexible operation for different products, increased worker safety, increased productivity, and reduced hazardous waste disposal costs.

ROSENMUND

Rosenmund Inc.
P.O. Box 668625
Charlotte, N.C. 28266-8625
704-398-1111 (Tel) 704-398-1089 (Fax)

Rosenmund Inc. serves the pharmaceutical, chemical, specialty chemical, and process industries with a full range of filters, filter-dryers, and Apovac[®] vacuum equipment. Engineered systems include solvent recovery, sterile and CIP designs, recirculation drying, heat transfer, reactor to filter-dryer packages, and customized applications.


ROSENMUND INC.

Filtration, Drying, Vacuum Equipment & Engineered Systems

INSTALLATION & EQUIPMENT HIGHLIGHTS


5M² Rosenmund Filter-Dryer[®]
Center Discharge-heated
paddle system

For details, please see reverse side.



Technical Specifications

5M² Rosenmund Filter-Dryer® Center Discharge-heated paddle system



Size	5M ²
Type	Paddle Filter-Dryer
Filtration Area	54 ft. ²
Working Slurry Volume	1530 gallons
Solids Volume	88 ft. ³
Design Pressure	50 psig/full vacuum
Design Temperature	-20° to 350°F
Material of Construction	Hastelloy C-276
	ASME Code
Agitator	Three-arm paddle, heated
Solids Discharge	Center
Filter Plate	Simple heated
	150 psig, 350°F
Jacket	Half-pipe coil
	150 psig, 350°F
Temperature Measurement	Thermowell
Dust Filter	Blowback, bag-type
Motors	40 HP, hydraulic, agitator rotation
	5 HP, hydraulic, all raising/lowering operations
Weight (empty)	26,000 pounds
Installation Date	August, 1990

Application Details

This pharmaceutical customer installed a 5M² Rosenmund Liquid Discharge Filter constructed in 316L stainless steel and a 5M² Rosenmund Filter-Dryer constructed in Hastelloy C-276. The first filter is employed for the filtration and washing of the two initial intermediate compounds which are both discharged as a slurry. The custom designed filter-dryer processes the finished product through a series of filtration, washing and drying operations. This installation, selected as an alternative to centrifuges and dryers, eliminates product handling, improves product yield, increases worker safety, and allows for sequential processing. Rosenmund supplied an APOVAC Vacuum/Gas Recirculation skid package, heat transfer module (HTM) and Kammerer Knife Gate Valve. A Rosenmund control system, based upon the Texas Instrument TI 305, operates each of the filters.

ROSENMUND

Rosenmund Inc.
P.O. Box 668625
Charlotte, N.C. 28266-8625
704-398-1111 (Tel) 704-398-1089 (Fax)

Rosenmund Inc. serves the pharmaceutical, chemical, specialty chemical, and process industries with a full range of filters, filter-dryers, and Apovac® vacuum equipment. Engineered systems include solvent recovery, sterile and CIP designs, recirculation drying, heat transfer, reactor to filter-dryer packages, and customized applications.


ROSENMUND INC.


Filtration, Drying, Vacuum Equipment & Engineered Systems

INSTALLATION & EQUIPMENT HIGHLIGHTS

8M² Rosenmund Filter-Dryer®
Center Discharge-heated
paddle system

For details, please see reverse side.





Size	8M ²
Type	Paddle Filter-Dryer
Filtration Area	86 ft. ²
Working Slurry Volume	2400 gallons
Solids Volume	153 ft. ³
Design Pressure	50 psig/full vacuum
Design Temperature	300°F
Material of Construction	316L stainless steel ASME Code
Agitator	Three-arm paddle, heated
Solids Discharge	Center
Filter Plate	Simple heated 90 psig, 300°F
Jacket	Half-pipe coil 90 psig, 300°F
Temperature Measurement ...	Thermowell
Dust Filter	Blowback, bag-type
Motors	40 HP, hydraulic, agitator rotation 5 HP, hydraulic, all raising/lowering operations
Weight (empty)	46,500 pounds
Installation Date	August, 1989

Application Details

The filter-dryer replaced existing centrifuge to dryer combinations at this pharmaceutical plant. The operator conducts filtration, smoothing, several displacement washes using different solvents and then dries this pharmaceutical intermediate to less than 0.1 percent moisture. Pre-drying is accomplished with initial hot gas recirculation using an Apovac® Vacuum/Gas Compressor System. Once the cake reaches its desired temperature vacuum drying begins with suction from above the cake through the dust filter and from below the cake through the filter media. The heated paddles stir the cake adding approximately 40 percent additional energy into the cake, which significantly shortens the drying cycle. The dried product is discharged into bins.

ROSENMUND

Rosenmund Inc.
P.O. Box 668625
Charlotte, N.C. 28266-8625
704-398-1111 (Tel) 704-398-1089 (Fax)

Rosenmund Inc. serves the pharmaceutical, chemical, specialty chemical, and process industries with a full range of filters, filter-dryers, and Apovac® vacuum equipment. Engineered systems include solvent recovery, sterile and CIP designs, recirculation drying, heat transfer, reactor to filter-dryer packages, and customized applications.

ROSENMUND INC.

Filtration, Drying, Vacuum Equipment & Engineered Systems

INSTALLATION & EQUIPMENT HIGHLIGHTS


10M² Rosenmund Filter-Dryer®
Center Discharge-heated
paddle system

For details, please see reverse side.



Technical Specifications

10M² Rosenmund Filter-Dryer[®] Center Discharge-heated paddle system



Size	10M ²
Type	Paddle Filter-Dryer
Filtration Area	108 ft. ²
Working Slurry Volume	3200 gallons
Solids Volume	180 ft. ³
Design Pressure	50 psig/full vacuum
Design Temperature	300°F
Material of Construction	316L stainless steel ASME Code
Agitator	Three-arm paddle, heated
Solids Discharge	Center
Filter Plate	Simple heated 90 psig, 300°F
Jacket	Half-pipe coil 90 psig, 300°F
Temperature Measurement ...	Thermowell
Dust Filter	Blowback, bag-type
Motors	40 HP, hydraulic, agitator rotation 5 HP, hydraulic, all raising/lowering operations
Weight (empty)	50,000 pounds
Installation Date	May, 1990

Application Details

This pharmaceutical installation conducts filtration, smoothing, washing and then dries this pharmaceutical intermediate to less than 1.0 percent moisture. To accomplish drying, the operator combines initial hot gas recirculation, vacuum drying through the dust filter and through the filtrate outlets, and stirring with the heated paddles. The heated paddles add approximately 40 percent additional energy into the cake, which significantly shortens the drying cycle. Rosenmund has completed the installation by supplying an APOVAC[®] Vacuum/Gas Compressor System, heat transfer module (HTM), fully automated PLC Control System, and a Kammerer Knife Gate Valve.

ROSENMUND

Rosenmund Inc.
P.O. Box 668625
Charlotte, N.C. 28266-8625
704-398-1111 (Tel) 704-398-1089 (Fax)

Rosenmund Inc. serves the pharmaceutical, chemical, specialty chemical, and process industries with a full range of filters, filter-dryers, and Apovac[®] vacuum equipment. Engineered systems include solvent recovery, sterile and CIP designs, recirculation drying, heat transfer, reactor to filter-dryer packages, and customized applications.

Combined Pressure Filter-Dryer Process Unit Improves Product Quality and Reduces Production Time

by Barry A. Perlmutter

A major U.S. pharmaceutical manufacturer produces two products through a series of reactions, filtration, washing, and drying steps in 12 and 18 hours, respectively. Current operation requires significant product handling as solids are transferred from floor to floor and building to building among reactors, centrifuges, and dryers. Bins, bags and carts are employed for movement. Test work demonstrated the feasibility of employing a single filter-dryer unit to replace existing equipment for reduced production time, improved yield and product quality, and minimized environmental and worker exposure.

This article begins with an introduction to Nutsche filtration and its use for the isolation of solid particles from liquid slurries. Topics include filtration pressure, elastic versus inelastic cake, cake heights, and filtration mechanisms. The article continues by examining product drying within a Nutsche filter. Drying techniques of vacuum or contact, convection, blow drying with hot gas, and gas recirculation are described. Environmental control strategies for volatile organic compounds and particulate dust are also discussed.

Test program methodology and test data from laboratory scale and pilot plant scale tests at the customer's site are illustrated. Data is included for filtration rates, washing cycles, and various drying techniques and sequences. Heat transfer coefficients are also calculated. The data is segmented for each product and then used to provide the necessary scale-up information for production models. The test work demonstrated the feasibility of employing a single filter-dryer unit to replace existing equipment for reduced production time, improved product quality, and minimized environmental and worker exposure. The paper concludes with a methodology for formulating a test program to analyze existing filtration, washing, and drying systems. Further, it shows how the data can be used to optimize operating variables, equipment and component selection, and overall system and process design.

Nutsche Filtration

Present day chemical processing requires a variety of solid-liquid separation techniques. One type of separation equipment is the classic agitated pressure Nutsche filter. This unit has been developed over the past 25 years from a simple pressure/vacuum filter into a sophisticated piece of equipment which is capable of performing many process operations including filtration, displacement washing, reslurry, smoothing and compression, drying, and product discharge either wet, dry or as a liquid.

Mechanisms of Filtration

Suspended solids are separated from a fluid stream by three mechanisms: inertial impaction, diffusional interception, and direct interception. For separating particles from a liquid, the predominate method is direct interception. In this case, the

particles are removed because they are larger than the pores or holes in the filter medium. In deep bed or cake filtration, the solid material itself becomes the filter medium. Small particles are removed by the tortuous path the fluid takes through the filter cake. These particles are removed as follows:

- In the real world most suspended particles, even if quite small when viewed from some directions, are irregular in shape and hence can "bridge" an opening.
- A bridging effect can also occur if two or more particles strike an opening simultaneously.
- Once a particle has been stopped by a pore, that pore is at least partially occluded, and subsequently will be able to separate even smaller particles from the liquid stream.

Filtration Pressure

The optimal pressure applied depends very much upon the specific particle characteristics of the product. Particles are generally amorphous, and vary in shape from spherical to irregularly shaped. These parameters will determine the type of cake formation and hence filtration pressure. Figure 1 shows filtration speed versus filtration pressure. Experience has shown that the majority of filtrations are in the range of 20 - 35 psig. Along this range, the filtration speed flattens for an inelastic




Figure 1. Filtration Rate Versus Filtration Pressure.

(incompressible) cake. The drop in filtration speed is even more dramatic for a compressible cake.

Cake Height

The feasible use of pressure filters begins with cake heights of approximately eight cm. In the majority of applications cake heights of 15 cm to 35 cm are used. Cake heights of larger than 50 cm are unusual, and filter sizing is always based upon cake height.

Filtration Media

The selection of the mesh size depends on the distribution of the particle sizes of the product and on the smallest particles present in a product, respectively.

The finest filter cloths have mesh sizes of around 20 microns. Considering the fact that the same filter cake has a certain filter effect, one accepts a product loss at the beginning of the filter cycle by using a larger, more rigid filter cloth (30 to 50 microns) which definitely offers a larger servicing period.

Alternatives to filter cloths include metal screens and multi-layer filter media. Metal filter cloths are more expensive than those made of woven polymer, but also have a corresponding longer service life. Multi-layer filter media are considered to be the solution to many filter problems. These are the most expensive but also provide a defined structure for long filter life, easy cleaning, and less vulnerability to damage.

Cake Washing

The ideal washing procedure is established by trials. Less porous, crystalline particles usually need a displacement washing while porous, lamina-like products are preferably washed through a reslurry washing. In many cases however a combination of both shows the best results.

Filter-Dryer Technology

In conventional processes the pressure filter is one link in a chain of processes.

Given that the filter, as a pressure vessel, is already of robust construction and is equipped with an agitator and quite often with heated walls, it is a logical progression to combine filtration with the following process step, drying. Product handling between operations is thus eliminated, and product loss, product contamination and environmental contamination are avoided.

Vacuum Contact Drying

To convert a pressure filter to a filter-dryer, the following modifications are necessary. First, heat transfer surfaces are examined. These include vessel walls, dome, base, filter plate

and agitator. The agitator system is especially critical as between 30 - 50% of the heat transfer occurs through the heated paddles. Where a "pasty" phase is expected during drying, the agitator construction and the drive unit are to be correspondingly strengthened.

The agitator paddles, in a filter-dryer, are designed as streamlined blades with a large surface area, each equipped with mixing/smoothing shovel or paddles.

Ideally the filter-dryer is fitted with a variable speed drive unit of either frequency controlled electric motor or hydraulic type. Should a viscous phase occur during drying, the vertically adjustable paddle offers the special feature of self-regulation of the cut into the cake to the point where the maximum admissible torque is reached. As a result, the agitator automatically tends to work into the "pasty" product, until the viscous stage is completed and the agitator is once more in the lowest position.

Next, the filter medium should be selected. During the drying stage, the agitator paddles move with relatively little clearance from the filter plate (five - 20 mm). This has the effect of placing more stress on the filter medium than is the case during simple filtration. For this reason, laminated, sintered woven metal wire filter plates (multimesh) are recommended. The higher cost of this type of filter plate is compensated by increased service life.

Convection Drying

In convection drying, the heat is primarily conveyed by a hot gaseous medium such as air or nitrogen. Generally, two methods of convection drying are employed. The first is blow drying.

After filtration, the cake is blown through from above with warm air or warm nitrogen. The humidity in the cake is absorbed and carried away by the warm gas. In this case the filter mesh serves also as a "dust filter." This method is only suitable, within certain limitations, for drying porous product-masses. With compressible or fine grained products, the gas flow, even under high pressure, soon ceases and the drying operation stops.

The second method is via hot gas recirculation. Gas recirculation is generally from above the cake (downflow). The cake can be agitated or remain static depending upon the drying characteristics.

Environmental Controls

Finally, there are several environmental concerns when one is considering a filter-dryer. These are product particles and volatile organic compound (VOC) solvent emissions. To capture solid particles, a vapor dust filter is incorporated. To recover

1. Rosenmund Filter-Dryer with connection for heating/cooling media, nozzles for suspension (a), wash-liquor (b), pressure gas (c), filtrate (d), product discharge (e), gas circuit (f)
2. Dust filter with heating jacket and filter hose cleaned with purge gas (g)
3. Pre-condensor
4. Ring liquid compressor
5. Ring liquid cooler
6. Ring liquid container with condensate or ring liquid overflow (h)
7. Exhaust gas cooler/condensator with exhaust (i) and gas circuit connection (f)
8. Gas heater (for gas circuit)




Figure 2. Schematic of the Rosenmund® Filter-Dryer Installation.

*A total of 15 trials were conducted for the
two pharmaceutical intermediates.*

the VOC solvent emissions, a system consisting of a liquid-ring vacuum pump (LRVP), a liquid-ring cooler, and an exhaust gas cooler, has proved to be very successful. As a rule the vaporized solvent is used as the liquid-ring. Should a high vacuum be required in order to achieve low moisture content, the vacuum pump can be preceded by a gas ejector.

In the case of the requirement for hot gas recirculation, a closed circuit with integrated condensor, blower and gas heater is employed. Here again the LRVP system has proved to be reliable, with the liquid-ring pump being used as compressor (with an approximate two bar pressure differential). The cooled gas stream from the condensor is re-heated to drying temperature in a gas heating unit. Figure 2 shows a typical layout.

The similarity of the condensation and compression system makes the described variations easy to incorporate. The capital cost of the combined installation is only slightly higher, and is limited to the additional bypass valves and piping required to switch over from vacuum to gas circuit and the reverse. Combined filter-dryer installations are particularly suitable for multi-purpose plants where products with diverse properties (temperature sensitivity, solvent coefficients, etc.) are to be processed.

Filter-Dryer Test Program

The pharmaceutical manufacturer has, as its objective in this test program, to eliminate the solids handling of their two pharmaceutical intermediates. Bins, bags and carts presently transfer the products from floor-to-floor and building-to-building, between reactors, centrifuges, and dryers. Meeting this objective will allow the operator to reduce production time from 12 and 18 hours, respectively, improve yield and product

quality, and minimize environmental and worker exposure.

Overview

On-site test work to eliminate solids handling was conducted using a Rosenmund® 0.6M² Center Discharge Paddle pilot filter and an Apovac® LRVP vacuum/gas compressor system. The test program provides data for optimizing filtration and washing times, cake height, drying times and heat transfer coefficients. The conclusion of the test program provides a guaranteed total cycle time as well as defined filter-dryer size and drying sequence.

Results

A total of fifteen trials were conducted for the two pharmaceutical intermediates. Table 1 details the type of data collected rather than the actual results themselves. The specific important results are described below. The tests gave an average of 0.03 - 0.09% LOD at a cake height for Product A of 42 centimeters and 45 centimeters for Product B. Bulk densities were calculated based upon dry solids weight just after initial filtration. Filtration and washing times were optimized based upon product quality concerns and cycle times.

Drying tests were done with heated jacket and base, using the Apovac® 136 to generate vacuum. Paddles were agitated but not heated. Before starting vacuum drying or agitation, the cake was blown with cold nitrogen for 15 minutes, followed by hot nitrogen recirculation with the Apovac for 30 minutes to heat the cake.

The calculation of the heat transfer coefficient "U" is shown in Table 2. Heat load is based on the quantity of water evaporated during vacuum drying. Time is based on vacuum drying only. Useable wall areas is based on measured cake depth. Product temperature is estimated to be halfway between the heating jacket temperature and the equilibrium temperature of water at the operating pressure. From the heat load, vessel wall and base area, and temperature difference between jacket and product, U is calculated. For the production model, U will be increased as heated paddle agitators will be incorporated along

A. Filtration

1. Filtrate, liters
2. Time, minutes
3. Filtration Rate, L/min-M²

B. Washing

1. Water Wash, liters
2. Time, minutes
3. Wash Rate, L/min-M²
4. Solvent Wash, liters
5. Time, minutes
6. Wash Rate, L/min-M²

C. Drying

1. Vacuum Drying Time, hours
2. Jacket Temperature, °F
3. Pressure, inches Hg absolute
4. Initial LOD, %
5. Final LOD, %
6. Yield, lbs. (dry)

D. Cake Depth, cm

Dry Solids, lbs.

Initial LOD %

Final LOD %

Liquid Evaporated, lbs.

Heat Load, Btu x 10⁵

Vacuum Drying Time, hours

Useable Wall Area, ft.²

Base Area, ft.²

Total Area, ft.²

Jacket Temperature, °F

Product Temperature, °F

Heat Transfer Coefficient: U = Btu/hr-ft.²-°F, is calculated from the above data.

Table 1. Pilot Test Data Collected on the 0.6M² Rosenmund® Center Discharge Paddle Filter-Dryer.

Table 2. Data Used to Calculate Heat Transfer Coefficients from Pilot Tests.

with the heated walls and base. The heated paddles will add approximately 25 percent to the total heat transfer area.

Total Liquid to Filter, liters
Liquid Pumping Rate L/minute
Time, minute
hour
Filtration/Washing Time, hour
Nitrogen Blow Through Time, hour
Drying Time, hour
Discharge Time, hour
Total Time, hours, is equal to the sum of the above times.

Table 3. Calculation of Batch Cycle Times.

Finally, Table 3 illustrates the calculation of the total cycle time for the process. It includes both slurry filling time and dry product discharge time. The single filter-dryer unit reduced the processing time of the two products by 22 and 38 percent, respectively.

Test Program Methodology

Whether one is considering a filter or filter-dryer, a well-defined test program is required. The following are the important items to be covered prior to undertaking such tests.

1. Set objectives for process.
2. Determine application and process designs.
3. Product and slurry information including particle size distribution, existing data, new product data, other particle characteristics, etc.
4. Cake thickness and filtration rates in order to maximize thickness and minimize filtration rates. This should include selection of the proper filter media.
5. Washing steps including displacement wash, reslurry wash, solvent types. The goal is to maximize product purity while minimizing wash liquid.

6. Drying requirements including final moisture content, type of drying, sequence of drying, agitated versus static cake, drying cycles, heat transfer coefficients, and dry product yield.
7. Available utilities, space, product, slurry, personnel, etc.
8. Unique requirements: cleanliness, toxic or hazardous, equipment modifications, material compatibility, insurance, other criteria.

Conclusions


The technology of Nutsche filtration has progressed from a simple un-agitated open filter to its current level of a filter-dryer. The integration of this type of equipment into a production process requires an innovative approach to problem-solving and a willingness to carefully examine existing process steps. This pharmaceutical manufacturer conducted a specific test program to meet predetermined process objectives. Filtration, displacement and reslurry, and drying rates and sequences were optimized in relation to product quality concerns. The result is the installation of a cost-effective agitated pressure filter-dryer which reduces production time, minimizes product handling, improves product quality and yield, and decreases environmental and worker exposures to organic solvent vapors and solid particulates.



Barry A. Perlmutter is National Sales and Marketing Manager of Rosenmund, Inc., responsible for direct sales, marketing and new business development for Rosenmund products. Mr. Perlmutter has more than 10 years of plant and process experience in the chemical, pharmaceutical, environmental and process industries. Previously, he was a senior marketing manager with Pall Corp. Mr. Perlmutter began his career as an Environmental Scientist with the U.S. Environmental Protection Agency. He has published and lectured extensively on the science of filtration. Mr. Perlmutter holds a BS in chemistry, MS in environmental engineering and an MBA. He is a member of ISPE and PDA.


Standard Filterplate

Non Heating




Heated Filterplate

Shell and tube
„Complex system“



Heated Filterplate

Non perforated filterplate
„Simple system“



ROSENMUND INC


Draw

6-9-88

0012

FILTRATION AND

DRYING SYSTEMS



TECHNOLOGY UPDATE

ROSENMUND INC.

Vol. 1, No. 3

September 1989

APOVAC® Vacuum Systems for Solvent Recovery Operations

The pharmaceutical and chemical process industries are concerned with controlling vapor emissions for several reasons: to meet environmental regulations, for increased worker safety, and to minimize product losses. Various methods can be used including flaring and carbon adsorption. These techniques are end-of-the-pipe solutions and are not integral with the process operation itself.

APOVAC Vacuum System


An improved technique for the recovery of solvents during vacuum distillation, evaporation, and drying operations is the APOVAC vacuum/gas compressor system. The skid-mounted engineered package is based upon a single stage liquid ring pump as shown in Figure 1. The system allows for the recovery and reuse of clean solvent with minimum air emissions and zero liquid emissions.

APOVAC Operation and Benefits

Condensation of the solvent vapors is carried out in the pump itself. The condensate and remaining vapors pass through a ring liquid cooler, specially designed to achieve maximum contact and cooling. The chilled liquid is collected in the ring

liquid tank for recirculation to the pump. The remaining vapors travel through an exhaust gas cooler where more vapors condense and flow counter-current back into the ring liquid tank. The exiting gas meets applicable solvent emission regulations. The process solvent, which is also the ring liquid, is recovered and reused with no liquid waste streams.

The vertical ring liquid cooler in the APOVAC design produces the lowest possible ring liquid temperature and therefore the lowest possible suction pressure. The design also cools the inert gas as it is forced through the ring liquid cooler tubes resulting in very low exit gas temperatures. These low gas temperatures, in conjunction with the exhaust gas cooler, provide for significantly greater condensation and solvent recovery than conventional vacuum recovery techniques. Finally, the APOVAC unit, because of its design, does not require an additional seal fluid circulation pump and control loop or interstage cooling, thus eliminating possible sources of leaks. The system, in conclusion, is low-maintenance and inherently safe to operate.



APOVAC 156

APOVAC System Applications

Rosenmund Filter® and Filter-Dryer®

The APOVAC unit is installed with the Rosenmund Filter for vacuum filtration and vacuum drying. The APOVAC is also employed for gas

recirculation drying and once-through gas blowing.

Drying, Distillation, Evaporation, Mixing and Reaction

The APOVAC unit is uniquely suited for these applications where

single or multiple solvents are used. Typical solvents include methylene chloride, isopropanol, acetone, hexane, toluene, propylene glycol and water.

Neutralization

Corrosive, toxic, and odorous vapors can be neutralized with the APOVAC unit. In this case, the ring liquid is a pH-adjusted caustic soda or similar liquid. A static mixer is also employed upstream of the liquid ring pump for high scrubbing efficiency and to protect the APOVAC from fast corrosion.*

Replacement of Existing Vacuum Systems

The APOVAC unit can replace existing steam jets where there is high water usage, environmental problems, or inefficiencies. For conventional oil vacuum pumps, replacement with an APOVAC unit can eliminate oil contamination problems as well as the high maintenance costs associated with these systems.

APOVAC Test Program

APOVAC units are available in several different designs depending upon the specific application and plant conditions. Units are also available for lease at your plant site. Challenge Rosenmund to solve your vacuum system and solvent emission problems today.




Figure 1: Schematic Representation of the APOVAC System

*For further details on the static mixer design, please see "Better Absorption, Try a Static Mixer", Rader, Mutsakis, Koch Engineering Co., Grosz-Ruell, Maugweiler, Sulzer Bros., Ltd., Chemical Engineering, July, 1989, p. 137-142.

Next Issue: We Want to Hear Your Ideas.


The next issue of TECH UPDATE will be published in December, 1989 and will focus on the Rosenmund procedure for conducting filtration and drying tests.

In addition, we are looking for engineers who would like to publish an article in TECH UPDATE. The article could relate to Rosenmund technology or to solving a specific production problem with a unique approach. We want to hear your ideas. Look for details in the December issue or use the Response Card to submit your ideas.

We look forward to seeing you at the Chem Show in New York City, December 4th-7th at our booth, number 1457.

ROSENMUND

Rosenmund, Inc.
P.O. Box 668625
Charlotte, North Carolina 28266-8625
704-398-1111



**To evacuate and
handle gases
and vapours**

**APOVAC®
Vacuum Systems**

with closed
ring-liquid circuit


SULZER BURCKHARDT

APOVAC® ANTIPOLLUTION VACUUM

APOVAC®

For the evacuation of gases and vapours with recovery of the resulting non-contaminated and reusable condensate, for neutralization of acid gases without polluting air and effluent.

Suction capacity up to 725 m³/h (415 cfm) for vacuum up to 40 mbar (28.7" Hg) with ejector up to 6 mbar (29.85" Hg).





Technical data

Type	Maximum suction volume without gas ejector		Ring liquid circulation		Coolant at $\Delta t = 4^\circ\text{C}/7^\circ\text{F}$		Standard motors IP54		Motors for compressor duty up to 2 bar abs		Nominal speed rpm	Weight with pump and motor approx.	
	m³/h	cfm	m³/h	gpm	m³/h	gpm	kW	HP	kW	HP		kg	lbs
136	84	50	0.25	1.1	1.0	4.4	4.0	7.5	4.5	7.5	3600	165	365
156	135	80	0.50	2.2	1.0	4.4	5.5	7.5	7	10	3600	190	420
2507	180	105	1.0	4.4	2.0	8.8	7.5	10	10	15	1800	380	830
2510	300	170	1.1	4.8	2.0	8.8	10	15	13	20	1800	430	950
3810	410	240	2.5	11	5.0	22	15	20	20	30	1200	1065	2350
3813	570	335	2.7	12	7.0	31	20	30	24	30	1200	1065	2350
3817	725	425	3.2	14	7.0	31	25	40	33	40	1200	1130	2500



Ask for dimension sheets and installation diagrams for the APOVAC system selected.

Vacuum pump types according to data in separate pump leaflets.



Standard equipment

- 21 Ring-liquid tank
- 22 Ring-liquid cooler
- 23 Liquid-ring vacuum pump
- 24 Discharge duct
- 25 Ring-liquid pipe with Flow Switch
- 26 Manually operated vacuum control valve
- 27 Non-return valve



Accessories

- 30 Ejector with or without bypass
- 31 Exhaust gas condenser
- 32 Flow indicator with switching unit
- 33 Level monitoring (or control)
- 34 Level sight glass
- 35 Ring-liquid thermometer

5084 5006

Table of connections

All flange connections ANSI 150 lbs RF 316 LSS	APOVAC 136 + 156	APOVAC 2507 + 2510	APOVAC 3810 + 3813	APOVAC 3817
(A) Suction branch	2"	2½"	3"	3"
(B) Gas discharge	2"	3"	4"	5"
(C) Filling connection	1"	1"	1"	1"
(D) Overflow	1"	1"	1"	1"
(E) Drain	1"	1"	1"	1"
(F) Coolant in/outlet	1"	1½"	2"	2"
(G) Coolant in/outlet	1"	1½"	2"	2"
(K) Mechanical seal leakage	½"	½"	½"	½"
(L) Level monitoring	1½" gas	1½" gas	1½" gas	1½" gas
(N) Temperature indication	1" gas	1" gas	1" gas	1" gas

Construction materials

— Standard version, all stainless steel, 316L.

Subject to dimensional and design alterations

T. Sobieszek

Environmentally-compatible Evacuation of Gases and Vapours

Reprint from the Sulzer Technical Review 4/1986
© 1987 by Sulzer Brothers Limited, Winterthur

2148

Liquid-ring pumps are suitable for handling all kinds of gases and vapours, and may be employed as vacuum pumps or compressors.

Liquid-ring pumps (*Fig. 1*) need no oil lubrication. They are rugged, simple and reliable. Once installed and put into operation, they can be virtually forgotten [1].

A detailed view of a mechanical pump assembly, likely a diaphragm pump. It features a grey cylindrical head, a red electric motor, and various piping and valves. The pump is mounted on a base and has several connection points for hoses or pipes.

1 Liquid-ring pump with IEC standard motor, non-return valve, separator and economy circuit.

With the help of a combination of liquid-ring pumps with gas ejectors and Roots blowers (COMBIVAC system) and by combining liquid-ring pumps with vapour ejectors, it is possible to operate in the medium-high vacuum range (*Figs. 5 and 6*).






Figure 10 is a line graph showing the relationship between Intake rate (m³/h) on the y-axis and Intake pressure (mbar_{abs}) on the x-axis. The y-axis ranges from 0 to 600 in increments of 100. The x-axis ranges from 0 to 1000 in increments of 200. The graph is divided into two regions: 'with gas ejector' (left, lower pressures) and 'without gas ejector' (right, higher pressures). There are eight sets of curves, each representing a different gas ejector type, labeled with their respective numbers: 3817, 3813, 3810, 2510, 2507, 156, 136, and 122. The curves for 'with gas ejector' are black, and the curves for 'without gas ejector' are pink. The curves show that the intake rate increases with intake pressure and is generally higher for ejectors with higher numbers. The 'with gas ejector' curves show a sharp increase in intake rate at low pressures, while the 'without gas ejector' curves show a more gradual increase at higher pressures.

SULZER TECHNICAL REVIEW 4/1986




4 Liquid-ring pump with upstream compact gas ejector and integrated automatic bypass valve.



5 Performance curves with working pressure ranges of standard COMBIVAC vacuum units.

Properties

The use of water as ring liquid in open systems is without problems, provided the water is not contaminated with corrosive, toxic and malodorous gases and vapours. If it is, the polluted water cannot be discharged to drainage but must be treated. In many cases water is not allowed as ring liquid for various reasons. Situations of this kind lead to a system with a closed ring liquid circulation (Fig. 7). On account of its environmental compatibility, we call this system APOVAC (an abbreviation from "Anti-Pollution Vacuum system"). As ring liquid this system can use practically any fluid occurring from a process in the form of condensate (Table 1). It is particularly ideal for recovering solvent vapours. Corrosive, toxic and malodorous gases and vapours can be neutralized. The cooling water does




7 APOVAC system.

not have to be treated again, since it is not polluted. The additional costs of such a system are soon offset by the economic benefits.

The variety of problems encountered in chemical process engineering calls for individual approaches, leading to answers best suited to the particular purpose. Just one example here is the

L-Vac system (Figs. 8 and 9), which does not have a downstream exhaust gas condenser, though of course it could be fitted with one. This system is of simpler layout than APOVAC. Both of them can accommodate specific user needs in the chemical industry (e.g. explosion-proof pump drives, corrosion-proof material for the pump, condenser and circulation tank or thermal insulation). The L-Vac unit is suitable for duties where recovery is not of prime importance, such as neutralizing acid gases with caustic soda solutions as ring liquid. Where the medium to be evacuated permits, it is also possible to operate with a water ring. As a rule the L-Vac is employed with ring liquids of high boiling point, and APOVAC with low-boiling ones.




6 Standard COMBIVAC Type 1000/S/2510, with integrated APOVAC system (intake capacity 800 m³/h, working pressure range 0.4 to 10 mbar, final vacuum 0.05 mbar).

Environmentally-compatible operation

The use of vacuum units with integrated closed ring liquid circuits is indicated wherever liquids or solvents are employed as product processing aids but are removed again before the intermediate or end product is finished. Recovery is desired wherever possible.

To ensure gentle product treatment, this is accomplished at temperatures



8 L-Vac system.

Table 2: Exhaust gas load as percentage of 1 kg air

Ring liquid	(°C)	Circulation without gas cooling	L-Vac	APOVAC
Acetone	-15	19.8	14.4	6.3
Methyl alcohol	-5	8.6	6.2	2.5
Toluene	+10	12.2	7.2	3.3
Butyl alcohol	+30	8.5	4.7	2.7

Table 1: Ring liquids used most in closed APOVAC circuits

Alcohols

Methyl alcohol	CH_3O
Ethyl alcohol	$\text{C}_2\text{H}_5\text{O}$
Propyl alcohol	$\text{C}_3\text{H}_7\text{O}$
Butyl alcohol	$\text{C}_4\text{H}_9\text{O}$
Octyl alcohol	$\text{C}_8\text{H}_{17}\text{O}$
Butyl alcohol	$\text{C}_4\text{H}_{10}\text{O}$
Ethylene glycol	$\text{C}_2\text{H}_4\text{O}_2$
Glycerine	$\text{C}_3\text{H}_8\text{O}_3$
Cyclohexanol	$\text{C}_6\text{H}_{12}\text{O}$
Benzyl alcohol	$\text{C}_7\text{H}_8\text{O}$

Phenols

Phenol	$\text{C}_6\text{H}_6\text{O}$
--------	--------------------------------

Ketones

Acetone	$\text{C}_3\text{H}_6\text{O}$
Methyl ethyl ketone	$\text{C}_4\text{H}_8\text{O}$

Ethers

Dimethyl ether	$\text{C}_2\text{H}_6\text{O}$
Diethyl ether	$\text{C}_4\text{H}_{10}\text{O}$
Furan	$\text{C}_4\text{H}_4\text{O}$

Inorganic compounds

Water	H_2O
Sodium hydroxide (caustic soda)	NaOH

Halogenated hydrocarbons

Methylene chloride	CH_2Cl_2
Chloroform	CHCl_3
1,1-dichloroethane	$\text{C}_2\text{H}_4\text{Cl}_2$
Trichloroethylene	C_2HCl_3
Chlorobenzene	$\text{C}_6\text{H}_5\text{Cl}$

Organic sulphur compounds

Methyl mercaptan	CH_3S
Ethyl mercaptan	$\text{C}_2\text{H}_5\text{S}$

Aromatics

Benzene	C_6H_6
Toluene	C_7H_8
Ethyl benzene	C_8H_{10}
Styrene (vinylbenzene)	C_8H_8

Others

High temperature oils	
Polyethylene glycol	

iliary fluids according to their own boiling characteristics.

From the aspect of recovery, it is desirable to provide a special vacuum unit for each process unit (mixer, dryer, filter, distillation column, reactor etc.) (Fig. 10) which:

- corresponds to the pressure ranges of the process,
- allows optimal recovery at the point of origin,
- causes the lowest possible exhaust gas load, in order to maximize recovery and minimize exhaust post-treatment costs.

If we examine the question of exhaust gas load due to solvents in the circulation systems described, assuming similar operating conditions, differences emerge which, today more than ever, must be taken into account when choosing a system (Table 2).

The conclusion to be drawn from this is that systems with efficient ring liquid and gas cooling cause 40 to 73% less exhaust gas load, allow higher recovery and accordingly demand less post-treatment outlay.

The immediate future will show how far future legislation to combat atmospheric pollution is going to affect the individual vacuum units at each process stage. The admissible ppm loads already fixed cannot be attained with conventional condensers alone by reason of the gas laws, and any further additional measures are ultimately a question of cost.

Ω



9 L-Vac Type 13/80 PMZF 2510 with partial ring liquid feed in pump suction branch for neutralization purposes.



10 APOVAC units installed in a separate room of a fine chemicals plant.

Bibliography

KÓCWIN, M.J. Liquid-Ring Pumps for Handling Gases and Vapours. Sulzer Technical Review (1973), No. 2, p. 63.

lower than those tolerated by the product itself, in vacuum ranges allowing evaporation and withdrawal of the aux-

SULZER BURCKHARDT

Sulzer-Burckhardt
Engineering Works Ltd
Dornacherstrasse 210
CH-4002 Basel
Telephone 061-35 00 20
Telegrams Motor Basel
Telex 62785 mbbs ch
Telefax 061-35 14 70

Sulzer-Burckhardt
Engineering Works Ltd
Winterthur Division
CH-8401 Winterthur
Telephone 052-81 11 22
Telex 896 060 62 szch
Telefax 052-23 15 44

Sulzer-Burckhardt

Sulzer Group

Product range

Activity program

Reciprocating process
compressors, lubricated or
non-lubricated, vertical or
balanced-opposed
Air-cooled high-pressure
compressors
Water-lubricated compressors
Liquid-ring vacuum pumps
Labyrinth-piston compressors

Hot-water and steam generating
plants
Components for nuclear
power stations
Gas turbines
Diesel engines
Marine equipment
Water turbines
Storage pumps and pump-turbines
Valves and control systems
for power stations
Syngas coolers
Pumps, compressors and blowers
Machinery for the fiber-cement
industry
Plants for the manufacture of paper
and millboard
Textile machines
Heating, ventilating and air conditioning
systems
Plumbing and fire-protection systems
Refrigerating and cryogenic plants,
heat pumps
Plants for distillation, evaporation
and drying
Centrifuges and mixing processes
Chemical and biotechnical
process installations
Water and wastewater treatment
plants
Water cooling systems
Irradiation plants
Medical engineering products
Locomotives and railcars
Electronic regulating and
control devices
Industrial gears
Parts manufacture
Foundry products
Engineering of complete plants




EMIL KAMMERER



KNIFE GATE VALVES from ROSENMUND



THE KAMMERER VALVE



SIZE	1"	2"	2½"	3"	4"	6"	8"	10"	12"	16"	20"	24"	32"	40"
L	11.02	16.81	19.76	22.28	25.94	33.86	42.72	49.80	61.42	80.31	90.55	111.81	146.85	O/R
L₁	2.87	4.88	5.95	6.97	7.99	10.94	14.17	18.11	20.98	28.74	33.46	39.13	51.97	O/R
L₂	5.71	10.24	12.32	14.25	16.46	22.36	28.86	36.18	42.44	56.69	66.0	81.14	107.87	O/R
H	3.15	3.94	4.72	4.72	4.72	4.72	5.91	5.91	5.91	7.87	7.87	7.87	9.84	O/R
D	4.53	6.5	7.28	7.87	8.66	11.22	13.39	16.73	19.09	22.83	26.38	30.7	40.35	O/R
K	3.12	4.75	5.5	6.0	7.5	9.5	11.75	14.25	17	21.25	25.0	29.5	36	O/R



28



THE KAMMERER VALVE



INTRODUCTION

The Kammerer Knife Gate Valve range has been developed to meet the wide variety of needs presented by today's industry for a totally enclosed valve to handle particulate solids in all forms and sizes.

DESCRIPTION

The Kammerer valve comprises a totally enclosed body containing a moveable slide plate (knife) into which is cut an orifice of equal size to the valve bore. In the open position this orifice is in line with the valve bore and presents an uninterrupted opening. In the closed position the slide plate completely seals the valve bore to give total shut off.

Operation of the slide plate can be manual, pneumatic, hydraulic or electro-mechanical.

APPLICATION

Kammerer valves are specifically designed for systems containing particulate solids. The range of sizes that can be handled extends from dusts to greater than $\frac{1}{4}$ ".


The simplest design is dust-tight across the slide plate while the most sophisticated design will hold vacuum to 10^{-5} torr or pressure to 600 psig.

SELECTION

Kammerer valves are divided into 5 main categories, each of which is described in detail. These are in increasing order of sophistication.

'G'
'GP'
'P'
'S'
'F'

These are described in the next two pages.




TYPE 'G'

Designed for dust-tight closure of process lines containing particulate solids where full bore opening and minimum valve thickness are required.

As in all Kammerer valves, the Type 'G' incorporates a hard steel slide plate with full bore orifice. The slide plate is located above a machined, hard steel or hard bronze guide ring. Mounted in the guide ring and in the upper housing above the knife are scraper rings to clean the slide plate. These scraper rings encircle the valve bore thus ensuring contact with the slide plate at all times. Close tolerance between plate and rings ensures that the plate is always clean. The hard steel slide plate guarantees a long service life. The Type 'G' valve is dust-tight across the slide plate but will withstand line pressures to 100 psig without leaking to the atmosphere.

The Type 'G' valve is recommended for fine, free-flowing solids. The valve should only be closed when the line is free of solids but will retain and prevent solids flow in the closed position.

The Type 'G' valve is available in sizes from 2" to 20" bore. As with all Kammerer valves, it can be closed in less than 2 seconds whether manually operated or fitted with an automatic operator.




TYPE 'GP'

For granular abrasive materials, the type 'GP' is preferred. This valve is also suitable for low pressure systems involving gases, aggressive vapors or liquids.

The Type 'GP', like the Type 'G', is pressure tight in the housing and stem opening to 100 psig while across the slide plate it will seal against 10 psig or 25 mbars vacuum absolute.


The Type 'GP' is especially suited for conveying systems where the slide plate has to cut through particles on closing. For this an O-ring seal is incorporated in the upper housing. Below this seal in the guide ring a small notch is cut. This allows granular material to pass the O-ring without damaging the elastomer (see below). Tests have shown that even after 2000 closures, the O-ring showed minimum abrasion.

Like the Type 'G', the 'GP' can be closed in less than 2 seconds even when manually operated.





THE KAMMERER VALVE



TYPE 'S' TYPE 'P'

For systems requiring a valve to seal against a pressure differential greater than 10 psig, then the Types 'P' and Type 'S' are recommended. Both valves share the same pressurized sealing system; only the method of applying pressure to the seal varies.


These valves are designed for systems where a very good seal is required against pressure or high vacuum. The Type 'P' and Type 'S' will seal vacuum greater than 1 torr absolute. Pressure differentials up to 60 psig across the slide plate can be sealed and the housing is tight to 150 psig.

Sealing is by means of a U-shaped elastomer seal which is pressurized once the slide plate is in the closed position and has ceased moving. The external pressure forces the U-seal against the stationary slide plate to create the required seal. At all other times the U-seal is relaxed away from the slide plate and there is therefore no wear of the U-seal due to the movement of the slide plate.

The Type 'P' is available manually operated only and incorporates a self contained system for pressurizing the U-seal.

The Type 'S' uses an external source such as compressed air to pressurize the U-seal.

Size details for the Type 'S' and Type 'P' are given on page 2.



TYPE 'F'

Where elastomers cannot be used, due to solvent attack or where the temperature exceeds 120°C (250°F) then the Type 'F' is preferred. The basic operation is identical to the Type 'S' or 'P' but the U-seal is replaced with a metallic piston located in the upper housing. The piston is sealed to the housing by means of O-rings and to the slide plate by means of Teflon seals. A solid Teflon piston can also be provided.

In the case of temperatures above 180°C (356°F), the Teflon can be replaced by high temperature materials. See Type 'FBK' on page 6.


The Type 'F', Type 'S' and Type 'P' share a common body configuration and are available in sizes from 1" to 40" bore as detailed on page 2. As with the Type 'S', the Type 'F' uses an external source such as compressed air to pressurize the piston seal.

The Type 'F' is particularly recommended for use with vacuum dryers where there is a constant exposure to solvent vapors at elevated temperatures.

STANDARD VALVE VARIATIONS

The preceding 5 types of Kammerer valves are standard items in the size ranges indicated. The normal material of construction is 'Silumin', a chemically inert alloy of Al, Si, and Mg. It combines light weight with excellent chemical and abrasion resistance. All valves are standard with high chrome stainless steel slide plates and stems.

If stainless steel is required for contact surfaces, then the valve bore can be lined as shown below. Hastelloy, Monel, Inconel, etc. can also be used.



Valves are also available fabricated totally from any alloy desired. (see page 6.)


'KII' AND 'KIII' VARIATIONS

The standard valves are generally designed to seal process lines in which there is little or no solids flow when the valve is being closed. They are also not normally designed for sealing lines containing liquids where the liquid is the carrier phase, although they will handle wet solids.

To meet these special needs, there are two variations on the standard design, the 'KII' and the 'KIII', described on page 5. Both variations are available with the 5 types of sealing systems ('G', 'GP', 'P', 'S' and 'F') and should normally be specified with both designations — i.e., 'GKII', 'SKIII', etc.



THE KAMMERER VALVE




TYPE 'KII'

When it is necessary to close the valve through a stationary column of solids or when solids are passing through the valve as it is closing, then the 'KII' designed is preferred.

This valve is available with all five forms of sealing systems ('G', 'GP', 'P', 'S', and 'F'). The valve body is modified to incorporate an angled slot which intersects the path of the closing slide plate. Any material carried into the valve body, as the slide plate closes, has the possibility to drop, via the slot, back into the downstream pipeline. Any material not recovered in this way is caught in a chamber which is an integral part of the valve housing. This chamber can be cleaned out periodically. The valve is also equipped with purge connections so that an external gas source can be introduced into the body cavities to blow them clear, while the slide plate is moving.

This valve is normally mounted in vertical lines but can also be used in some inclined locations. It is not recommended for horizontal lines.

The 'KII' valve is available in all the same sizes as the standard valves, although dimensions may vary slightly. Materials of construction are the same as the standard valves.



TYPE 'KIII'

When handling a liquid, especially a viscous liquid, the standard Kammerer valves are not normally suitable as the liquid can enter into, and remain in, all the body cavities as the slide plate opens or closes. For this reason the 'KIII' design was developed to eliminate all cavities where liquids could be trapped.


Essentially, the 'KIII' valve comprises the upper part of a standard Kammerer valve body, with a tapered hopper section welded below. The slide plate and spindle are exposed on the lower side to the hopper section. A cast, machined center section locates the slide plate and contains the sealing mechanism which can again be Type 'G', 'GP', 'P', 'S' or 'F'. The hopper section is straight-sided on its two long sides and tapered 45° at the ends. All crevices are eliminated and the valve can be used on the most viscous of liquids without leaving any product in the valve.

The 'KIII' valve is only mounted in vertical lines. Sizes range from 1" NB up to 40" NB. Materials of construction are carbon steel, stainless steel, alloys, etc.

SPECIAL VALVES

Kammerer valves are also available in many special designs to meet requirements not covered by the valves mentioned previously. As there are many variations, any special requirements should be referred to ROSENMUND INC. so that the best design can be submitted.

Two special valve designs have however become 'standardized' and are described below.



TYPE 'HD'


Specifically for high pressure use from 120 psig to 600 psig pressure differential across the slide plate, the Type 'HD' valve modifies the concept of the Type 'F'. The standard square piston is replaced with a thinner rectangular piston with a profiled surface indenting a gasket material located in the slide plate. The slimmer design of the piston allows greater pressures to be applied against the gasket, thus ensuring a good seal at high pressure.

Two gaskets are normally mounted in the slide plate; one around the orifice to seal the valve body when the valve is open, and one in the solid portion of the slide plate to seal the valve against the differential pressure when the valve is in the closed position.

The Type 'HD' is available in all standard sizes and can be fabricated in all weldable materials.



THE KAMMERER VALVE



TYPE 'FBK'

The Type 'FBK' valve is designed specifically for high temperature applications above 500°F where there is also a need to seal against a pressure differential greater than 15 psig. The valve utilizes the same basic design as the Type 'F' but also includes aspects of the Type 'HD'. Temperatures up to 1500°F and pressures up to 450 psig have been successfully handled using this design.

All elastomeric and plastics materials are eliminated. The piston incorporates a profiled lower surface which indents into a high temperature gasket material set into the slide plate, as in the Type 'HD'. The pressure gas moving the piston is sealed from the valve body by means of a flexible metallic bellows welded to the upper surface of the bellows.

Again the 'FBK' design is available in all standard sizes and in all common weldable materials (subject to the temperature requirements).

Your Local Representative:

MATERIALS OF CONSTRUCTION

Body	Silumin (Anodized surface finish also available) Cast Iron Cast Steel Stainless Steel Hastelloy Other Weldable Alloys
Slide Plate	High Chrome Stainless Steel Hastelloy Other Alloys
Stem	Stainless Steel Hastelloy Other Alloys
Guide Ring	Chemical Hard Bronze Stainless Steel Other Alloys
Seals	EDPM Buna N® Silicon Rubber Viton® Teflon® Grafoil® Italon®
Valve Bore Lining	Any Metal

RANGE

The following represent the limits obtained with Kammerer valves to date:

Pressure.....to 600 psig

Vacuum.....to 10⁻⁵ torr abs

Temperature.....to 1500°F

We specialize in making valves for individual requirements. Kindly indicate your requirements and we will be pleased to build your valve.

SIZE

Standard sizes are featured on page 2. For sizes not featured, please contact ROSENMUND INC.

Valves are available flanged ANSI 150, ANSI 300, BS, DIN 2502 or DIN 2632. Please contact ROSENMUND INC. for full details.

OTHER PRODUCTS

Also available in the ROSENMUND range of Chemical Process Equipment:

- The ROSENMUND® Filter
- NUTREX® Reactor Filter Dryer
- KUEHNI Liquid Extraction System

For further details on all KAMMERER products, contact

ROSENMUND INC.

ROSENMUND INC.

CHEMICAL PROCESS EQUIPMENT

2969 INTERSTATE STREET

P. O. BOX 668625

CHARLOTTE, NC 28266-8625

TELEPHONE (704) 398-1111

TELEX 802079 Rosenmund cha

CHEMICAL PROCESS EQUIPMENT

CHARLOTTE, NC 28208 USA

802079 ROSENMUND CHA